

anti- ARMETL1 antibody

Product Information

Catalog No.: FNab00593

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: arginine-rich, mutated in early stage tumors-like 1

Synonyms: None
Observed MW: 21 kDa
UniprotID: Q49AH0

Application

Reactivity: Human

Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:200-1:1000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab00593(ARMETL1 Antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00593(ARMETL1 Antibody) at dilution of 1:300