

anti- BCL11A antibody

Product Information

Catalog No.:	FNab00837
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Functions as a myeloid and B-cell proto-oncogene. May play important roles in leukemogenesis and hematopoiesis. An essential factor in lymphopoiesis, is required for B-cell formation in fetal liver. May function as a modulator of the transcriptional repression activity of ARP1(By similarity).

Immunogen information

Immunogen:	B-cell CLL/lymphoma 11A(zinc finger protein)	
Synonyms:	CTIP1, EVI9, KIAA1809, ZNF856	
Observed MW:	None	
Uniprot ID :	Q9H165	

Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, IHC	
Recommended dilution: IHC: 1:20-1:200		


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889



Image:



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab00837(BCL11A antibody) at dilution of 1:100

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com