

anti- CD11B/Integrin alpha M antibody

Product Information

Catalog No.:	FNab01411
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	integrin, alpha M(complement component 3 receptor 3 subunit)
Synonyms:	CD11B, CR 3 alpha chain, CR3A, Integrin alpha M, ITGAM, MAC 1, MAC1A, MO1A, Neutrophil adherence receptor, SLEB6
Observed MW:	None
Uniprot ID :	P11215

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01411 (CD11B/Integrin alpha M antibody) at dilution of 1:800