
anti- CD43 antibody

Product Information

Catalog No.:	FNab01481
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6G11
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	sialophorin
Synonyms:	GALGP,SPN
Observed MW:	115 kDa
UniprotID :	P16150

Application


Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01481(SPN Antibody) at dilution of 1:100 . Heat mediated antigen retrieved with Tris-EDTA buffer, pH9.0



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01481(SPN Antibody) at dilution of 1:1000