

anti- Collagen Type II antibody

Product Information

Catalog No.:	FNab01837
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	collagen, type II, alpha 1
Synonyms:	Alpha 1 type II collagen, ANFH, AOM, COL11A3, COL2A1, Collagen alpha 1(II) chain, Collagen Type II, collagen, type II, alpha 1, SEDC
Observed MW:	141 kDa
UniprotID :	P02458

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human lung tissue slide using FNab01837(COL2A1 Antibody) at dilution of 1:50



Mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab01837(COL2A1 antibody) at dilution of 1:500