
anti- DNAJB3 antibody

Product Information

Catalog No.:	FNab02449
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May operate as a co-chaperone of the male germ cell-and haploid stage-specific Hsp70 proteins.

Immunogen information

Immunogen:	DnaJ(Hsp40) homolog, subfamily B, member 3
Synonyms:	HCG3
Observed MW:	28kd
UniprotID :	Q8WWF6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:1000-1:4000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis using FNab02449(DNAJB3 antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02449(DNAJB3 antibody) at dilution of 1:800