

## anti- EEF1B2 antibody

### Product Information

Catalog No.:	FNab02645
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8A11
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

EF-1-beta and EF-1-delta stimulate the exchange of GDP bound to EF-1-alpha to GTP.

### Immunogen information

Immunogen:	eukaryotic translation elongation factor 1 beta 2
Synonyms:	EEF1B, EF1B
Observed MW:	32 kDa
Uniprot ID :	P24534

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:4000; IF: 1:50-1:500
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02645( EEF1B2 Antibody) at dilution of 1:2000