

## anti- FXR1 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab03256  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | fragile X mental retardation, autosomal homolog 1 |
| Synonyms:    | FXR1, hFXR1p                                      |
| Observed MW: | 60-80 kDa   |
| Uniprot ID : | P51114  |

### Application


|                       |   |
|-----------------------|---|
| Reactivity:           | Human, Mouse, Rat   |
| Tested Application:   | ELISA, WB, IHC, IP, IF  |
| Recommended dilution: | WB: 1:1000-1:10000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100 |
| Image:                |   |

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human kidney using FNab03256(FXR1 antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab03256(FXR1 antibody) at dilution of 1:300