

anti- FXR2 antibody

Product Information

Catalog No.:	FNab03257
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

RNA-binding protein.

Immunogen information

Immunogen:	fragile X mental retardation, autosomal homolog 2
Synonyms:	FMR1L2
Observed MW:	95 kDa
UniprotID :	P51116

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-FXR2 (IP:FNab03257, 4ug;
Detection:FNab03257 1:500) with HeLa cells
lysate 3200ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab03257(FXR2 antibody)
at dilution of 1:500