

## anti- HJURP antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab03892   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres. Acts as a specific chaperone for CENPA and is required for the incorporation of newly synthesized CENPA molecules into nucleosomes at replicated centromeres. Prevents CENPA-H4 tetramerization and prevents premature DNA binding by the CENPA-H4 tetramer. Directly binds Holliday junctions.

### Immunogen information

|              |                                       |
|--------------|---------------------------------------|
| Immunogen:   | Holliday junction recognition protein |
| Synonyms:    | FAKTS, FLEG1, URLC9                   |
| Observed MW: | 70kd,83-85kd                          |
| UniprotID :  | Q8NCD3                                |

### Application

|                          |              |
|--------------------------|--------------|
| Reactivity:              | Human        |
| Tested Application:      | ELISA, WB    |
| Recommended dilution:WB: | 1:500-1:2000 |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Image:



A549 cells were subjected to SDS PAGE followed by western blot with FNab03892(HJURP antibody) at dilution of 1:800

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)