

anti- L-Plastin antibody

Product Information

Catalog No.:	FNab04833
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	lymphocyte cytosolic protein 1(L-plastin)
Synonyms:	None
Observed MW:	68 kDa
UniprotID :	P13796

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:6000; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

A431 cells were subjected to SDS PAGE followed by western blot with FNab04833(L-Plastin antibody) at dilution of 1:3000

