

## anti- LSM1 antibody

### Product Information

Catalog No.:	FNab04870
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Plays a role in replication-dependent histone mRNA degradation. Binds specifically to the 3'-terminal U-tract of U6 snRNA.

### Immunogen information

Immunogen:	LSM1 homolog, U6 small nuclear RNA associated(S. cerevisiae)
Synonyms:	CASM
Observed MW:	18-20 kDa
UniprotID :	O15116

### Application


Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:3000; IHC: 1:50 - 1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab04870(LSM1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04870(LSM1 antibody) at dilution of 1:600