

anti- MIA antibody

Product Information

Catalog No.: FNab05173

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

Immunogen information

Immunogen: melanoma inhibitory activity

Synonyms: None
Observed MW: None
Uniprot ID: Q16674

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human skin cancer slide using FNab05173(MIA Antibody) at dilution of 1:50