

anti- MYL4 antibody

Product Information

Catalog No.:	FNab05488
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Regulatory light chain of myosin. Does not bind calcium.

Immunogen information

Immunogen:	myosin, light chain 4, alkali; atrial, embryonic
Synonyms:	ALC1, AMLC, GT1, MLC1, MYL4, Myosin light chain 4, PRO1957
Observed MW:	22kd
Uniprot ID :	P12829

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



human liver tissue were subjected to SDS PAGE followed by western blot with FNab05488(MYL4 antibody) at dilution of 1:1000