

anti- ORAI3 antibody

Product Information

Catalog No.:	FNab06004
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Key regulator or component of store-operated Ca(2+) channel and transcription factor NFAT nuclear import.

Immunogen information

Immunogen:	ORAI calcium release-activated calcium modulator 3
Synonyms:	TMEM142C
Observed MW:	~30 kDa
UniprotID :	Q9BRQ5

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human endometrial cancer tissue slide using FNab06004(ORAI3 Antibody) at dilution of 1:50



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab06004(ORAI3 Antibody) at dilution of 1:300