

anti- PAX8 antibody

Product Information

Catalog No.:	FNab06175
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.

Immunogen information

Immunogen:	paired box 8
Synonyms:	paired box 8, Paired box protein Pax 8, PAX8
Observed MW:	48 kDa, 58 kDa
UniprotID :	Q06710

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:2000-1:10000; IHC: 1:1000-1:5000; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab06175(PAX8 Antibody) at dilution of 1:1200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab06175(PAX8 antibody) at dilution of 1:1000