

## anti- Pre-IL18 antibody

### Product Information

Catalog No.:	FNab06765
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4A7
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	interleukin 18(interferon-gamma-inducing factor)
Synonyms:	None
Observed MW:	24 kDa
Uniprot ID :	Q14116

### Application

Reactivity:	Human, Hamster
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:10-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)

HeLa cells were subjected to SDS PAGE followed by western blot with FNab06765( Pre-IL18 Antibody) at dilution of 1:1000

