

## anti- SAA1 antibody

### Product Information

Catalog No.:	FNab07572
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Major acute phase protein.

### Immunogen information

Immunogen:	serum amyloid A1
Synonyms:	PIG4, SAA, SAA1, Serum amyloid A protein, serum amyloid A1, TP53I4
Observed MW:	14 kDa
UniprotID :	P0DJ18

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemical of paraffin-embedded human liver using FNab07572(SAA1 antibody) at dilution of 1:50



Rat liver tissue were subjected to SDS PAGE followed by western blot with FNab07572(SAA1 antibody) at dilution of 1:1000