

## anti- TPPP3 antibody

### Product Information

Catalog No.:	FNab08900
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

### Immunogen information

Immunogen:	tubulin polymerization-promoting protein family member 3
Synonyms:	CGI-38, p20, p25gamma, TPPP/p20
Observed MW:	19 kDa
UniprotID :	Q9BW30

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human colon carcinoma using FNab08900( TPPP3 Antibody) at dilution of 1:100



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab08900( TPPP3 Antibody) at dilution of 1:1000