

## anti- TSPAN33 antibody

### Product Information

Catalog No.:	FNab09063
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Plays an important role in normal erythropoiesis. It has a role in the differentiation of erythroid progenitors(By similarity).

### Immunogen information

Immunogen:	tetraspanin 33
Synonyms:	HPen, PEN, Penumbra, Proerythroblast new membrane, tetraspanin 33, Tspan 33, TSPAN33
Observed MW:	32 kDa, Observed 70 kDa
UniprotID :	Q86UF1

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab09063(TSPAN33 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09063(TSPAN33 antibody) at dilution of 1:400