

anti- WDR46 antibody

Product Information

Catalog No.:	FNab09496
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Scaffold component of the nucleolar structure. Required for localization of DDX21 and NCL to the granular compartment of the nucleolus.

Immunogen information

Immunogen:	WD repeat domain 46
Synonyms:	BING4, C6orf11, FP221, WD repeat domain 46, WDR46
Observed MW:	75 kDa
Uniprot ID :	O15213

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IP: 1:200-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09496(WDR46 antibody) at dilution of 1:400