

anti- XKR6 antibody

Product Information

Catalog No.:	FNab09544
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	XK, Kell blood group complex subunit-related family, member 6
Synonyms:	XKR6, XRG6, XK related protein 6, C8orf7, C8orf5
Observed MW:	42 kDa
Uniprot ID :	Q5GH73

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09544(XKR6 Antibody) at dilution of 1:1000