

---

## anti- HBc antibody

### Product Information

Catalog No.:	FNab09775
Size:	1mg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
IsoType:	IgG1
Storage:	PBS with 0.1% Proclin300 pH 7.3,store at 2-8° C up to 21 days, and at -20° C for long-term storage. (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen: hepatitis B virus core antigen

### Application

Reactivity: Recommended for ELISA, RIA and CLIA. Each laboratory should determine an optimum working titer for use in its particular platform application.