

anti- STAG3 Antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab09779 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

None

Immunogen information

| | |
|--------------|---|
| Immunogen: | stromal antigen 3 |
| Synonyms: | Cohesin subunit SA 3, SCC3 homolog 3, STAG3, stromal antigen 3, Stromalin 3 |
| Observed MW: | 155 kDa |
| UniprotID : | Q9UJ98 |

Application


| | |
|-----------------------|---|
| Reactivity: | Human |
| Tested Application: | ELISA, IP, IHC, IF |
| Recommended dilution: | IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:50-1:500 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab09779(STAG3 Antibody) at dilution of 1:200



IP Result of anti-STAG3 (IP: FNab09779, 4ug; Detection: FNab09779 1:300) with human testis tissue lysate 4000ug.