
anti- T4 antibody

Product Information

Catalog No.:	FNab09809
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C7
IsoType:	IgG1
Storage:	20mMPBS, pH7.4 with 0.1% sodium azide, 2-8°C

Immunogen information

Immunogen:	T4
------------	----

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com