

## anti- CD74 antibody

### Product Information

Catalog No.:	FNab09829
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2C5
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	CD74
Synonyms:	CD74, DHLAG, HLADG, Ia GAMMA, Ii, p33
Observed MW:	34 kDa
Uniprot ID :	P04233

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:6000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Raji cells were subjected to SDS PAGE followed by western blot with FNab09829( CD74 Antibody) at dilution of 1:3000