

## anti- MTNR1A antibody

### Product Information

Catalog No.:	FNab10523
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Melatonin receptor type 1A
Synonyms:	Mel-1A-R, Mell1a receptor, Melatonin receptor type 1A
Observed MW:	39kd
Uniprot ID :	P48039

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded mouse brain slide using FNab10523(MTNR1A antibody) at dilution of 1:100



Rat kidney were subjected to SDS PAGE followed by western blot with FNab10523(FNab10523 antibody) at dilution of 1:1000