

# anti- ALDH1A1 antibody

## **Product Information**

Catalog No.:	FNab00284
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1G9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Binds free retinal and cellular retinol-binding protein-bound retinal. Can convert/oxidize retinaldehyde to retinoic acid(By similarity).

#### **Immunogen information**

Immunogen:	aldehyde dehydrogenase 1 family, member A1
Synonyms:	ALDC, ALDH E1, ALDH1, ALDH11, ALDH1A1, ALHDII, PUMB1, RALDH 1, RALDH1, Retinal dehydrogenase 1
Observed MW:	55 kDa
Uniprot ID :	P00352

# Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC, IF		
Recommended dilution: WB: 1:200-1:1000; IHC: 1:50-1:500; IF: 1:10-1:100			
Image:			

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan,<br/>Hubei, China(430206)Tel :( 0086)027-87384275Fax: (0086)027-87800889www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab00284(ALDH1A1 Antibody) at dilution of 1:200



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00284(ALDH1A1 antibody) at dilution of 1:500

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)