

# anti- APOL1 antibody

#### **Product Information**

Catalog No.: FNab00503

Size: 100μg Form: liquid

Purification: Protein A+G purification

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse

Clonality: monoclonal

Clone ID: 5E2
IsoType: IgG2a

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

May play a role in lipid exchange and transport throughout the body. May participate in reverse cholesterol transport from peripheral cells to the liver.

### **Immunogen information**

Immunogen: apolipoprotein L, 1

Synonyms: APO L, APOL, APOL I, APOL1, Apolipoprotein L, Apolipoprotein L I,

apolipoprotein L, 1, Apolipoprotein L1

Observed MW: 44 kDa UniprotID: 014791

## **Application**

Reactivity: Human

Tested Application: ELISA, WB, IHC, FC

Recommended dilution: WB: 1:10000-1:50000; IHC: 1:50-1:500

Image:


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1





Immunohistochemistry of paraffin-embedded human liver using FNab00503(APOL1 antibody) at dilution of 1:50



1 ug human plasma tissue was subjected to SDS PAGE followed by western blot with FNab00503(APOL1 Antibody) at dilution of 1:50000