

# anti- APOL4 antibody

## **Product Information**

Catalog No.:	FNab00505
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a role in lipid exchange and transport throughout the body. May participate in reverse cholesterol transport from peripheral cells to the liver(By similarity).

#### **Immunogen** information

Immunogen:	apolipoprotein L, 4
Synonyms:	APOL IV, APOL4, Apolipoprotein L IV, apolipoprotein L, 4, Apolipoprotein L4, APOLIV
Observed MW:	50-52 kDa
Uniprot ID :	Q9BPW4

## Application


Reactivity:	Mouse, Human	
Tested Application:	ELISA, IP, WB	
Recommended dilution: WB: 1:500-1:2000; IP: 1:200-1:1000		
Image:		

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





IP Result of anti-APOL4 (IP:FNab00505, 4ug; Detection:FNab00505 1:1000) with HepG2 cells lysate 4000ug.

HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00505(APOL4 Antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com