

# anti- ARAP3 antibody

#### **Product Information**

Catalog No.: FNab00523

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

# **Background**

Phosphatidylinositol 3,4,5-trisphosphate-dependent GTPase-activating protein that modulates actin cytoskeleton remodeling by regulating ARF and RHO family members. Is activated by phosphatidylinositol 3,4,5-trisphosphate(PtdIns(3,4,5)P3) binding. Can be activated by phosphatidylinositol 3,4-bisphosphate(PtdIns(3,4,5)P2) binding, albeit with lower efficiency. Acts on ARF6, RAC1, RHOA and CDC42. Plays a role in the internalization of anthrax toxin.

#### **Immunogen information**

Immunogen: ArfGAP with RhoGAP domain, ankyrin repeat and PH domain 3

Synonyms: CENTD3
Observed MW: 45 kDa

Uniprot ID: Q8WWN8

## **Application**

Reactivity: Human

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:20-1:200

## Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



# Image:



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab00523(ARAP3 Antibody) at dilution of 1:50

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>