

anti- ARC/ARG3.1 antibody

Product Information

Catalog No.: FNab00524

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: activity-regulated cytoskeleton-associated protein

Synonyms: None

Observed MW: 45-50 kDa UniprotID: Q7LC44

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC, IP

Recommended dilution: WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>






Immunohistochemistry of paraffin-embedded human brain using FNab00524(ARC antibody) at dilution of 1:100



IP Result of anti-ARC/ARG3.1 (IP:FNab00524, 4ug; Detection:FNab00524 1:300) with mouse brain tissue lysate 4000ug.





SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab00524(ARC antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com