

anti- ARD1A antibody

Product Information

Catalog No.: FNab00526

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: ARD1 homolog A, N-acetyltransferase(S. cerevisiae)

Synonyms: None
Observed MW: 26 kDa
UniprotID: Q6P4J0

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human breast cancer using FNab00526(ARD1A antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00526(ARD1A antibody) at dilution of 1:500