

# anti- ARF4 antibody

#### **Product Information**

Catalog No.: FNab00532

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking; may modulate vesicle budding and uncoating within the Golgi apparatus.

# Immunogen information

Immunogen: ADP-ribosylation factor 4

Synonyms: ARF2
Observed MW: 21 kDa
UniprotID: P18085

## **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC, IP, IF, FC

Recommended dilution: WB: 1:500-1:2000; IP: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1






Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab00532(ARF4 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



IP Result of anti-ARF4 (IP:FNab00532, 3ug; Detection:FNab00532 1:600) with mouse brain tissue lysate 5200ug.





mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00532(ARF4 Antibody) at dilution of 1:800

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com