

# anti- ARF5 antibody

#### **Product Information**

Catalog No.: FNab00533

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

# **Background**

GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking; may modulate vesicle budding and uncoating within the Golgi apparatus.

## Immunogen information

Immunogen: ADP-ribosylation factor 5

Synonyms: None
Observed MW: 21 kDa
UniprotID: P84085

# **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


### 1

## Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human breast cancer using FNab00533(ARF5 antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00533(ARF5 antibody) at dilution of 1:500