

# anti- ARHGAP15 antibody

#### **Product Information**

Catalog No.: FNab00547

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state. Has activity toward RAC1. Overexpression results in an increase in actin stress fibers and cell contraction.

### Immunogen information

Immunogen: Rho GTPase activating protein 15

Synonyms: ARHGAP15, BM 030, BM 046, BM046

Observed MW: 65-70kd Uniprot ID: Q53QZ3

# **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


# 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human lymphoma using FNab00547(ARHGAP15 antibody) at dilution of 1:50



human liver tissue were subjected to SDS PAGE followed by western blot with FNab00547(ARHGAP15 antibody) at dilution of 1:500