

# anti- ARPP-19 antibody

### **Product Information**

Catalog No.:	FNab00606
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

#### **Immunogen information**

Immunogen:	cyclic AMP phosphoprotein, 19 kD
Synonyms:	None
Observed MW:	19 kDa
Uniprot ID :	P56211


## Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, IHC, WB, IP, IF		
Recommended dilution: WB: 1:500-1:1000; IF: 1:10-1:100 ; IHC: 1:20-1:200			
Image:			

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





Immunohistochemistry of paraffin-embedded human gliomas using FNab00606(ARPP-19 antibody) at dilution of 1:50

IP Result of anti-ARPP-19 (IP:FNab00606, 3ug; Detection:FNab00606 1:500) with mouse brain tissue lysate 5200ug.




# 2

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com





Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00606(ARPP-19 antibody) at dilution of 1:500

3

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com