

anti- ARSA antibody

Product Information

Catalog No.:	FNab00610	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

Hydrolyzes cerebroside sulfate.

Immunogen information

Immunogen:	arylsulfatase A
Synonyms:	ARSA, arylsulfatase A, ASA, Cerebroside sulfatase, MLD
Observed MW:	56kd
Uniprot ID :	P15289

Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200			
Image:			

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com





Immunohistochemistry of paraffin-embedded human cervical cancer using FNab00610(ARSA antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00610(ARSA antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

2