

# anti- ARSK-Specific antibody

### **Product Information**

Catalog No.:	FNab00613
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

#### **Immunogen information**

Immunogen:	arylsulfatase family, member K
Synonyms:	None
Observed MW:	55 kDa
Uniprot ID :	Q6UWY0


## Application

Reactivity:	Human			
Tested Application:	ELISA, WB			
Recommended dilution:WB: 1:500-1:2000				
Image:				

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





HeLa cells were subjected to SDS PAGE followed by western blot with FNab00613(ARSK-Specific antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com