

# anti- ASL antibody

### **Product Information**

Catalog No.: FNab00641

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen: argininosuccinate lyase

Synonyms: argininosuccinate lyase, Arginosuccinase, ASAL, ASL

Observed MW: 49 kDa Uniprot ID: P04424

## **Application**

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:1000; IHC: 1:100-1:500

Image:


#### 1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab00641( ASL Antibody) at dilution of 1:200



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00641( ASL Antibody) at dilution of 1:500