

anti- ASNS antibody

Product Information

| Catalog No.: | FNab00644 |
|---------------|---|
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | \geq 95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |
| Background | |

Background

None

Immunogen information

| Immunogen: | asparagine synthetase |
|--------------|-----------------------------------|
| Synonyms: | ASNS, asparagine synthetase, TS11 |
| Observed MW: | 55-64 kDa |
| UniprotID : | P08243 |

Application


| Reactivity: | Human, Mouse, Rat | | |
|--|-------------------|--|--|
| Tested Application: | ELISA, WB, IP | | |
| Recommended dilution: WB: 1:500-1:2000; IP: 1:500-1:1000 | | | |
| Image: | | | |

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





IP Result of anti-ASNS (IP:FNab00644, 4ug; Detection:FNab00644 1:500) with HeLa cells lysate 1200ug.

HeLa cells were subjected to SDS PAGE followed by western blot with FNab00644(ASNS antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com