

## anti- ASS1 antibody

### Product Information

Catalog No.:	FNab00650
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4C3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Is indirectly involved in the control of blood pressure.

### Immunogen information

Immunogen:	argininosuccinatesynthetase 1
Synonyms:	Argininosuccinate synthase, argininosuccinatesynthetase 1, ASS, ASS1, Citrulline aspartate ligase, CTLN1
Observed MW:	45 kDa
UniprotID :	P00966

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:100-1:500; IF: 1:20-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human liver using FNab00650(ASS1 antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab00650(ASS1 Antibody) at dilution of 1:2000