

anti- ATGL antibody

Product Information

Catalog No.: FNab00680

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: patatin-like phospholipase domain containing 2

Synonyms: 1110001C14Rik, Adipose Triglyceride Lipase, ATGL, DESNUTRIN,

FP17548, IPLA2 zeta, PEDF R, PNPLA2, Transport secretion protein 2,

TTS 2.2, TTS2, TTS2.2

Observed MW: 55 kDa Uniprot ID: Q96AD5

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:1000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





A431 cells were subjected to SDS PAGE followed by western blot with FNab00680(ATGL antibody) at dilution of 1:500