

anti- ATL3 antibody

Product Information

Catalog No.: FNab00684

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

GTPase tethering membranes through formation of trans-homooligomers and mediating homotypic fusion of endoplasmic reticulum membranes. Functions in endoplasmic reticulum tubular network biogenesis.

Immunogen information

Immunogen: atlastinGTPase 3

Synonyms: ATL3, Atlastin 3, atlastinGTPase 3

Observed MW: 61 kDa
UniprotID: Q6DD88

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IF, WB

Recommended dilution: WB: 1:500-1:3000; IF: 1:10-1:100

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00684(ATL3 Antibody) at dilution of 1:1000