

anti- ATN1 antibody

Product Information

Catalog No.:	FNab00685
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcriptional corepressor. Recruits NR2E1 to repress transcription. Promotes vascular smooth cell(VSMC) migration and orientation(By similarity). Corepressor of MTG8 transcriptional repression. Has some intrinsic repression activity which is independent of the number of poly-Asn(polyQ) repeats.

Immunogen information

Immunogen:	atrophin 1
Synonyms:	D12S755E, DRPLA
Observed MW:	105-110 kDa
Uniprot ID :	P54259

Application


Reactivity:	Human, Mouse	
Tested Application:	ELISA, WB, IHC	
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200		

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com



Image:



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab00685(ATN1 Antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00685(ATN1 Antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>