

# anti- ATRIP antibody

### **Product Information**

Catalog No.:	FNab00738
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Required for checkpoint signaling after DNA damage.Required for ATR expression, possibly by stabilizing the protein.

## Immunogen information

Immunogen:	ATR interacting protein
Synonyms:	AGS1, ATR interacting protein, ATRIP
Observed MW:	86-90kd
UniprotID :	Q8WXE1

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:2000
Image:	

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





IP Result of anti-ATRIP (IP:FNab00738, 3ug; Detection:FNab00738 1:1000) with MCF-7 cells lysate 2800ug.

MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00738(ATRIP antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>