

# anti- ATXN2L antibody

#### **Product Information**

Catalog No.: FNab00741

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Involved in the regulation of stress granule and P-body formation.

## **Immunogen information**

Immunogen: ataxin 2-like

Synonyms: A2D, A2LG, A2LP, A2RP

Observed MW: 140 kDa
Uniprot ID: Q8WWM7

## **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00741(ATXN2L Antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00741(ATXN2L Antibody) at dilution of 1:1000