

anti- ATXN3L antibody

Product Information

Catalog No.:	FNab00743
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Deubiquitinating enzyme that cleaves both 'Lys-48'-linked and 'Lys-63'-linked poly-ubiquitin chains(in vitro).

Immunogen information

Immunogen:	ataxin 3-like
Synonyms:	ataxin 3 like, ATX3L, ATXN3L, MJDL, Putative ataxin 3 like protein
Observed MW:	41 kDa
Uniprot ID :	Q9H3M9

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00743(ATXN3L Antibody) at dilution of 1:500