

anti- Aurora A antibody

Product Information

Catalog No.: FNab00746

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: aurora kinase A

Synonyms: AIK, ARK 1, ARK1, AURA, AURKA, Aurora A, aurora kinase A,

Aurora related kinase 1, Aurora/IPL1 related kinase 1, AURORA2, Breast tumor amplified kinase, BTAK, hARK1, STK15, STK6, STK7

46 kDa

Observed MW: 46 kDa Uniprot ID: 014965

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:500-1:1000; IF: 1:50-1:200; IHC: 1:100-1:400

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab00746(AURKA Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00746(AURKA Antibody) at dilution of 1:600