

anti- B4GALT4 antibody

Product Information

Catalog No.: FNab00775

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.

Immunogen information

Immunogen: UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 4

Synonyms: B4Gal T4, B4GALT4, Beta 1,4 GalTase 4, beta4Gal T4

Observed MW: 50 kDa
Uniprot ID: 060513

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab00775(B4GALT4 antibody) at dilution of 1:500